

MARLIN® rapid

Qualities

- short absorption time
- braided
- synthetic
- coated
- ensures light knot sliding
- excellent tissue penetration
- made from polyglycolic acid (PGA)

Marlin®rapid degrades in the human body by hydrolysis with very low tissue reaction. Tensile strength is reduced by 50 per cent after 7 days, after about 42 days degradation is completed.

Application

- subcutaneous
- intracutaneous
- paediatric surgery
- gastro intestinal
- mucosa
- ligatures
- urological
- genital
- fascial
- peritoneum
- parenchymal
- bronchus and lung
- ophthalmic










How supplied

- needled sutures
- ligatures

MARLIN® rapid

synthetic, absorbable
braided, coated









Needled sutures

| Needle description | Article-No. | Colour | USP | EP | Length | UoS |
|---|-------------|--------|-----|-----|--------|-----|
|  HR 17 | 17261074 | undyed | 5/0 | 1 | 70 cm | 24 |
| | 17261082 | undyed | 4/0 | 1,5 | 70 cm | 24 |
| | 17261090 | undyed | 3/0 | 2 | 70 cm | 24 |
|  HR 22 | 17261140 | undyed | 4/0 | 1,5 | 70 cm | 24 |
| | 17261144 | undyed | 3/0 | 2 | 70 cm | 24 |
| | 17261150 | undyed | 2/0 | 3 | 70 cm | 24 |
|  HR 26 | 17261183 | undyed | 3/0 | 2 | 70 cm | 24 |
| | 17261186 | undyed | 2/0 | 3 | 70 cm | 24 |
| | 17261192 | undyed | 0 | 3,5 | 70 cm | 24 |
|  HR 27s | 17261257 | undyed | 2/0 | 3 | 70 cm | 24 |
| | 17261266 | undyed | 0 | 3,5 | 70 cm | 24 |
| | 17261272 | undyed | 1 | 4 | 70 cm | 24 |
| | 17261275 | undyed | 2 | 5 | 70 cm | 24 |
|  HR 37s | 17261412 | undyed | 2/0 | 3 | 90 cm | 24 |
| | 17261420 | undyed | 0 | 3,5 | 90 cm | 24 |
| | 17261422 | undyed | 1 | 4 | 90 cm | 24 |
| | 17261425 | undyed | 2 | 5 | 90 cm | 24 |
|  HR 43s | 17261511 | undyed | 0 | 3,5 | 90 cm | 24 |
| | 17261515 | undyed | 1 | 4 | 90 cm | 24 |
| | 17261521 | undyed | 2 | 5 | 90 cm | 24 |
|  HR 48s | 17261566 | undyed | 1 | 4 | 90 cm | 24 |
| | 17261572 | undyed | 2 | 5 | 90 cm | 24 |
|  HRT 26 | 17261676 | undyed | 3/0 | 2 | 45 cm | 24 |
|  HRT 37 | 17261700 | undyed | 2/0 | 3 | 90 cm | 24 |
| | 17261705 | undyed | 0 | 3,5 | 90 cm | 24 |
| | 17261710 | undyed | 1 | 4 | 90 cm | 24 |

MARLIN® rapid

synthetic, absorbable
braided, coated

Needled sutures

| Needle description | Article-No. | Colour | USP | EP | Length | UoS |
|---|--|--------------------------------------|--------------------------|------------------|----------------------------------|----------------------|
|  HS 15 | 17261844 17261848 | undyed undyed | 4/0 3/0 | 1,5 2 | 45 cm 45 cm | 24 24 |
|  HS 20 | 17261879 17261883 | undyed undyed | 3/0 2/0 | 2 3 | 45 cm 45 cm | 24 24 |
|  DS 19 | 17262030 17262034 17262049 | undyed undyed undyed | 4/0 3/0 3/0 | 1,5 2 2 | 45 cm 45 cm 70 cm | 24 24 24 |
|  DS 24 | 17262066 17262091 17262070 | undyed undyed undyed | 3/0 3/0 2/0 | 2 2 3 | 45 cm 70 cm 45 cm | 24 24 24 |
|  DSM 11 | 17262107 | undyed | 6/0 | 0,7 | 45 cm | 24 |
|  DSM 13 | 17262130 17262134 | undyed undyed | 5/0 4/0 | 1 1,5 | 45 cm 45 cm | 24 24 |
|  DSM 19 | 17262160 17262164 | undyed undyed | 4/0 3/0 | 1,5 2 | 45 cm 45 cm | 24 24 |
|  DSM 24 | 17262181 17262190 17262185 17262194 | undyed undyed undyed undyed | 3/0 3/0 2/0 2/0 | 2 2 3 3 | 45 cm 70 cm 45 cm 70 cm | 24 24 24 24 |

MARLIN® rapid

synthetic, absorbable
braided, coated

Ligatures

| | Article-No. | Colour | USP | EP | Length | UoS |
|-----------|-------------|--------|-----|-----|---------|-----|
| Ligatures | 17260017 | undyed | 3/0 | 2 | 3x45 cm | 24 |
| | 17260020 | undyed | 2/0 | 3 | 3x45 cm | 24 |
| | 17260022 | undyed | 0 | 3,5 | 3x45 cm | 24 |
| | 17260070 | undyed | 3/0 | 2 | 6x45 cm | 24 |
| | 17260074 | undyed | 2/0 | 3 | 6x45 cm | 24 |
| | 17260081 | undyed | 0 | 3,5 | 6x45 cm | 24 |
| | 17260085 | undyed | 1 | 4 | 6x45 cm | 24 |
| | 17260227 | undyed | 2/0 | 3 | 2x70 cm | 24 |
| | 17260230 | undyed | 0 | 3,5 | 2x70 cm | 24 |
| | 17260236 | undyed | 1 | 4 | 2x70 cm | 24 |
| | 17260240 | undyed | 2 | 5 | 2x70 cm | 24 |